

FIT TEST REPORT

10/19/2020

ID NUMBER GSF1
LAST NAME Schwind
FIRST NAME John
COMPANY Global Safety First, LLC
LOCATION Sea Girt NJ

TEST DATE 06/06/2020 12:24
DUE DATE 06/06/2021

PORTACOUNT S/N 8048201912
N95-COMPANION™ Y

RESPIRATOR Global Safety First, LLC N1901L Filtering Facepiec [100] N95
MANUFACTURER Global Safety First, LLC
MODEL N1901L
RESPIRATOR STYLE Filtering Facepiece
RESPIRATOR SIZE Large

PROTOCOL OSHA 29CFR1910.134
PASS LEVEL 100
APPROVAL
EFFICIENCY<99% Y

<u>EXERCISE</u>	<u>DURATION (sec)</u>	<u>FIT FACTOR</u>	<u>PASS</u>
NORMAL BREATHING	60	200+	Y
DEEP BREATHING	60	200+	Y
HEAD SIDE TO SIDE	60	200+	Y
HEAD UP AND DOWN	60	200+	Y
TALKING	60	200+	Y
GRIMACE	15	excl.	
BENDING OVER	60	200+	Y
NORMAL BREATHING	60	200+	Y
OVERALL FF		200+	Y

FITTEST OPERATOR _____
JAS

DATE _____

NAME _____
John Schwind

DATE _____

NOTE: