



ReadiMask Sizing Guidance

The Small size fits people with smaller faces and/or children. Many of our nurses are rather petite and they prefer the small size.

The dimensions are 6 and 7/8 inches wide x 5 and 3/8 inches high.

The Large size is really a medium/large and fits most adults both men and women.

The dimensions are 7 and 7/8 inches wide by 5 and 7/8 inches high.

The Extra Large size fits people with larger than average faces and features.

The dimensions are 8 and 7/8 inches wide by 7 inches high.

The height measurement is more critical than the width, so if you measure from the bridge of your nose to the bottom of your chin, then add 1 inch, that would be the minimum height you would need.